

UNIT COUNTING GUIDE




CLASPS
 Akers, 1 unit each
 RPI, RPA, RPR, and Equipoise
 1½ Units each



HORIZONTAL FINISH LINES
 3 to 6 teeth including retention
 1½ units, 7 or more teeth or acrylic
 palate, 2 units.




CLASPS WITH LUGS
 1½ Units each



KENNEDY BAR
 ¼ unit per tooth passed



RESTS
 Plain rest ¼ unit
 Individual res with strut ½ unit



APRON OR LINGUAL PLATE
 ¼ unit per tooth covered, plus lingual
 bar. Cingulum bar ½ unit per tooth




LINGUAL BAR
 For each tooth passed ¼ unit
 2 units maximum, 1 unit minimum



HORSESHOE
 For each tooth passed, ¼ unit
 2 units maximum



STRENGTHENING BARS
 1½ units
 Mesh reinforcement 2 units
 Double strengthening bar 2 units



CIRCULAR BAR
 Double palatal bars 2 units
 for each bar



DISTAL RETENTIONS
 With vertical finish line 1 unit




PALATAL BARS
 2 units, wide palatal strap 3 units



HORIZONTAL FINISH LINE
 Up to 2 teeth including retention
 1 unit



LOWER CAST BASE
 4 units



UPPER CAST BASE
 5½ units. With mesh post dam
 6½ units

Turn Around Times

Frames = Next Day

Frames with Acetal Resin Clasp = Next Day + 1

Frames with Metal Teeth, Tube Teeth, Onlays, or Bite Openers = 4 Days or less

Flexibles = 4 Days